



LOW TEMPERATURE WATER BATH

Environment Friendly Engineering Solution Company



In Association with SVCH-Technologii, Moscow (Russia)

www.kerone.com

ISO 9001-2008 | ISO 9001-2015 | EMS 14001 | OHSAS 18001

About Kerone

KERONE now renowned name in serving specialized need of customers with best quality and economical process Heating /cooling and drying products, manufactured in high quality environment by well trained and qualified workforce(special purpose machineries) .

KERONE is pioneer in application and implementation engineering with its vast experience and team of professionals. KERONE is devoted to serve the industry to optimize their operations both economically and environmentally with its specialized heating and drying solutions.


“Enhance the value of customer operation through our customer need centric engineering solution”

Mission

- To enhance the value of customer operation through our customer need centric engineering solution.
- We are committed to provide our customers, unique and best in class products in Industrial heating, drying and cooling segment, with strategic tie-up for the technical know-how with renowned leader in the industry specific segment.

Vision

- Turn into world leader in providing specialized, top-notch quality and ecological industrial heating, cooling and drying solution across the globe.
- To attain global recognition as best of quality and environment friendly engineering solution company.



Introduction to

Low Temperature Water Bath

A Low Temperature water bath is a device that retains water at a constant temperature. It is mainly used in laboratory for incubations at the lower temperature. KERONE designs and manufactures a Low Temperature Water Bath both standard and customized for the specialized need; every Low Temperature Water Bath manufactured by KERONE is conformance with the ISO standards for the application in Medical, Agricultural, Bio-Chemical Industrial, Research Laboratories & various other Industries.

KERONE offer fully controllable Low Temperature Water Bath which can be easily controlled and maintained desired temperature for the specified duration of the time.

[kerone/low-temperature-water-bath](http://kerone.com/low-temperature-water-bath)

www.kerone.com/low-temperature-water-bath.php



Advantages

- Dual LED screen displays temperature and time
- Protection against compressor overload, high/low voltage and current
- Offers precise temperature control and uniformity
- Inner and outer chamber made of stainless steel is corrosion resistant and sturdy
- Equipped with caster for easy installation
- Convenient and reliable operation
- Stainless steel working chamber
- Intelligent microcomputer temperature adjuster
- Digital display of temperature setup and measure
- PID parameters controlling function
- High accuracy and dependability of temperature controlling
- Cooling device assembly

3D View

Trusted Partner



Our Clients



KERONE

Kerone Engineering Solutions Limited

UNIT I

B/10, Marudhar Industrial Estate, Goddev Fatak road, Bhayander(E), Mumbai-401105

EMAIL

info@kerone.com | sales@kerone.com
marketing@kerone.com

UNIT II

Plot No. B-47, Addl. MIDC Anandnagar, Ambernath (East), Dist. Thane- 421506

WEBSITE

www.kerone.com | www.kerone.net
www.keroneindia.com

UNIT III

Plat NO. W-104, Addl. MIDC Anandnagar, Ambernath East, Thane- 421506. (India)

Phone

+91-22-28150612/13/14
+91-251-2620542/43/44/45/46