

# WATER STILL (DISTILLATION UNIT)

Environment Friendly Engineering Solution Company



In Association with SVCH-Technologii, Moscow (Russia)

[www.kerone.com](http://www.kerone.com)

ISO 9001-2008 | ISO 9001-2015 | EMS 14001 | OHSAS 18001

# About Kerone

KERONE now renowned name in serving specialized need of customers with best quality and economical process Heating /cooling and drying products, manufactured in high quality environment by well trained and qualified workforce(special purpose machineries) .

KERONE is pioneer in application and implementation engineering with its vast experience and team of professionals. KERONE is devoted to serve the industry to optimize their operations both economically and environmentally with its specialized heating and drying solutions.


**“Enhance the value of customer operation through our customer need centric engineering solution”**

## *Mission*

- To enhance the value of customer operation through our customer need centric engineering solution.
- We are committed to provide our customers, unique and best in class products in Industrial heating, drying and cooling segment, with strategic tie-up for the technical know-how with renowned leader in the industry specific segment.

## *Vision*

- Turn into world leader in providing specialized, top-notch quality and ecological industrial heating, cooling and drying solution across the globe.
- To attain global recognition as best of quality and environment friendly engineering solution company.



## Introduction to **Water Still (Distillation Unit)**

KERONE holds in proficiency in providing high quality, efficient and sustainable solutions with its range of laboratory equipments such as water still or water distillation unit used in various laboratories and industrial applications both standard and customized. The water still or water distillation units are precisely designed to meet the client's application requirements and comply with the standards. The water still or water distillation units are constructed with very high quality stainless still to avoid any corrosive reaction in both internal and exterior of body.

We also manufacture glass water distillation unit which is a unique & intelligent water distillation unit provides great user experience. The heating element is ceramic bobbin heater with Nichrome heating wire. The water still or water distillation unit has very low conductivity (less than 1 mho). We offer it in majorly in two variant models 1.5lt/hr distillation water output for continuous run & 4.5 liters/hr capacity as standard off the shelf product, including this we also design and build various custom build water still or water distillation units.

**[kerone/water-still-distillation-unit](http://kerone.com/water-still-distillation-unit)**

[www.kerone.com/water-still-distillation-unit.php](http://www.kerone.com/water-still-distillation-unit.php)



# Features

- Strong construction
- Rust proof
- Easy to install
- Mantle Type heating system
- Low maintenance
- Both 3 phase/single phase electrical power supply model
- Intelligent & user Friendly.

# Applications

The water still or water distillation units find its application in testing of:

- For Distillation of Water (General Lab Use)
- Steel Manufacturing Units
- Ferrous and Alloy Steel Foundries
- Forging Units
- Wire manufacturing Units
- QA and QC testing laboratories.



# 3D View

# Trusted Partner



# Our Clients





**KERONE**

Kerone Engineering Solutions Limited

### UNIT I

B/10, Marudhar Industrial Estate, Goddev Fatak road, Bhayander(E), Mumbai-401105

#### EMAIL

info@kerone.com | sales@kerone.com  
marketing@kerone.com

### UNIT II

Plot No. B-47, Addl. MIDC Anandnagar, Ambernath (East), Dist. Thane- 421506

#### WEBSITE

www.kerone.com | www.kerone.net  
www.keroneindia.com

### UNIT III

Plat NO. W-104, Addl. MIDC Anandnagar, Ambernath East, Thane- 421506. (India)

#### Phone

+91-22-28150612/13/14  
+91-251-2620542/43/44/45/46